

+  
WE  
KNOW  
ENDO.



# EndoActivator<sup>®</sup>

An easy to use & safe solution  
to activate your irrigation

---

 Dentsply  
Sirona  
Endodontics

# Safe. Effective. Simple.

Contra-angled design for easy access to posterior teeth

Light-touch ON/OFF button

Ergonomic, cordless, and battery-operated handpiece



Quick connect/disconnect Activators for various procedures

3-speed sonic motor


Designed by Clifford J. Ruddle DDS, Robert H. Sharp DDS, Pierre Machtou DDS

## One fact

Fluid activation, in well shaped canals, plays an important role in debridement and disinfection of the root canal system.

## One solution

The EndoActivator® is designed to safely and vigorously energize intracanal irrigants during endodontic treatment.



Courtesy of Clifford J. Ruddle DDS

## Multiple benefits

### Safe

- Strong, flexible medical grade polymer tips
- Single patient use
- Uncoated and non-cutting tips

### Effective

- Creates fluid hydrodynamics
- Improves debridement and the disruption of the smear layer and biofilm

### Simple

- Very simple clinical technique
- Intuitive device
- Ideal in practice when portability is required

## SYSTEM KIT INCLUDES

1x EndoActivator® Handpiece/Driver (Battery Included)	REF A 0911
Activator Tips:	
25x Small (#15/.02),	
25x Medium (#25/.04)	X
25x Large (#35/.04)	
100 x Protective Barriers	REF A 0914
1x Directions for Use	X
1x Sequence Card	7
1x Instructional DVD	X
Possibility to order each part separately	

## SPARE PARTS



**Small**  
#15/.02  
REF A 0913 022 015 00



**Medium**  
#25/.04  
REF A 0913 022 025 00



**Large**  
#35/.04  
REF A 0913 022 035 00