

+
**WE
KNOW
ENDO.**

MAILLEFER

Conform Fit™ Gutta-Percha

The right fit and flow you can feel

 **Dentsply
Sirona**



Next-generation gutta-percha. Made possible by Conform Fit™ technology.

Gutta-percha master cones should ideally fit the shapes created by your files. Yet, traditional hand-rolled gutta-percha points have made it hard to achieve a precise match. The next level of precision performance is here — gutta-percha master cones with Conform Fit technology.



The right fit you can feel. Variable tapered cones provide a snug fit that reaches the apex for the right tug back.



Improved flow to the end. Conform Fit gutta-percha has a better flow and transfers heat up to 5 mm beyond the heat source.



Enhanced patient safety. Master cones are natural rubber latex-free.



The simplicity of a total solution. Precision-matched shapes and sizes remove the guesswork.



Can be used at lower working temperature. Lower temperatures reduce the risk of harm to the periodontal ligament (PDL).




Convenient tab and colour markings for easy identification.




The right fit you can feel


You need a consistent and confident cone fit. However finding the right tug back with hand-rolled gutta-percha points can be tough. Traditional gutta-percha may seem to fit yet can actually be loose at the apex. Conform Fit™ gutta-percha cones feature a variable taper with a 40% more consistent shape than traditional gutta-percha. This creates a snug apical fit and an accurate tug back that you can feel.



Improved flow to the end

The unique gutta-percha formulation heats more consistently for better flow and allows deeper thermal penetration into the material. Heat transfers through the master cone up to 5 mm beyond the heat source—so you can be sure you have warm, flowable gutta-percha all the way to the apex. This enables the clinician to utilise a lower heat plugger temperature, thereby improving patient safety and achieving a complete fill.





Variable tapered cones provide a snug fit that reaches the apex for the right tug back.

Improved flow at a lower working temperature and transfers heat up to 5 mm beyond the heat source.*

40% more consistent shape and fit* than traditional hand-rolled gutta-percha points.

Master cones are natural rubber latex-free.

Part of a system-based solution for simplicity.

Convenient tab for easy identification.

*Source: Internal testing. Data available.




Enhanced patient safety

Conform Fit™ gutta-percha is not made with natural rubber latex. Our advanced formulation can be used at a lower working temperature, reducing the risk of harm to the periodontal ligament (PDL), which can be damaged by higher heat.



The simplicity of a total solution

Gutta-percha master cones with Conform Fit technology are part of a total system solution that takes the guesswork out of achieving a precise fit. Use them with corresponding files, gutta-percha points and paper points with confidence. The innovative tab design simplifies identification.





Part of a total system solution

WaveOne® Gold Conform Fit™ Gutta-Percha


Color	Size	Pack	Item #
Yellow	Small	60-pack	B00WGGPF00SML
Red	Primary	60-pack	B00WGGPF00PRI
Green	Medium	60-pack	B00WGGPF00MED
White	Large	60-pack	B00WGGPF00LRG
Assorted	S/P/M/L	60-pack	B00WGGPF00AST



WaveOne® Gold Reciprocating Files and WaveOne® Gold Paper Points

ProTaper Gold® Conform Fit™ Gutta-Percha


Color	Size	Pack	Item #
Yellow	F1	60-pack	B00PGGPF000F1
Red	F2	60-pack	B00PGGPF000F2
Blue	F3	60-pack	B00PGGPF000F3
Assorted	F2/F3	60-pack	B00PGGPF00F23
Assorted	F4/F5	60-pack	B00PGGPF00F45



ProTaper Gold® Endodontic Files and ProTaper Gold® Paper Points

ProTaper Next® Conform Fit™ Gutta-Percha


Color	Size	Pack	Item #
Red	X2	60-pack	B00PNGPF000X2
Blue	X3	60-pack	B00PNGPF000X3
Assorted	X2/X3	60-pack	B00PNGPF00X23
Assorted	X4/X5	60-pack	B00PNGPF00X45



ProTaper Next® Endodontic Files and ProTaper Next® Paper Points

TruNatomy™ Conform Fit™ Gutta-Percha

Color	Size	Pack	Item #
Yellow	Small	60-pack	B00TNGPF00SML
Red	Prime	60-pack	B00TNGPF00PRM
Green	Medium	60-pack	B00TNGPF00MED
Assorted	S/P/M	60-pack	B00TNGPF00AST



TruNatomy™ Endodontic Files and TruNatomy™ Paper Points

